*Flavobacterium columnare*

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- Bader, J. A. and E. B. Shotts, Jr. (1998): Identification of flavobacterium and flexibacter species by species-specific polymerase chain reaction primers to the 16s ribosomal rna gene. *J. Aquat. Anim. Health*, 10, 311-319.
- Bader, J. A., K. E. Nusbaum and C. A. Shoemaker (2003a): Comparative challenge model of *Flavobacterium columnare* using abraded and unabraded channel catfish, *Ictalurus punctatus* (Rafinesque). *J. Fish Dis.*, 26, 461-467.
- Bader, J. A., C. A. Shoemaker and P. H. Klesius (2003b): Rapid detection of columnaris disease in channel catfish (*Ictalurus punctatus*) with a new species-specific 16-S rRNA gene-based PCR primer for *Flavobacterium columnare*. *J. Microbiol. Methods*, 52, 209-220.
- Darwish, A. M., A. A. Ismaiel, J. C. Newton and J. Tang (2004): Identification of *Flavobacterium columnare* by a species-specific polymerase chain reaction and renaming of ATCC43622 strain to *Flavobacterium johnsoniae*. *Mol. Cell. Probes*, 18, 421-427.
- Panangala, V. S., C. A. Shoemaker and P. H. Klesius (2007): TaqMan real-time polymerase chain reaction assay for rapid detection of *Flavobacterium columnare*. *Aquacult. Res.*, 38, 508-517.
- Suebsing, R., R. Suebsing, J. Kampeera, J. Kampeera, S. Sirithammajak, S. Sirithammajak, B. Withyachumnarnkul, B. Withyachumnarnkul, W. Turner and W. Kiatpathomchai (2015): Colorimetric method of loop-mediated isothermal amplification with the pre-addition of calcein for detecting *flavobacterium columnare* and its assessment in tilapia farms. *J. Aquat. Anim. Health*, 27, 38-44.
- Toyama, T., K. Kita-Tsukamoto and H. Wakabayashi (1996): Identification of *Flexibacter maritimus*, *Flavobacterium branchiophilum* and *Cytophaga columnaris* by PCR targeted 16s ribosomal DNA. *Fish Pathol.*, 31, 25-31.
- Triyanto, A. Kumamaru and H. Wakabayashi (1999): The use of PCR targeted 16S rDNA for identification of genomovars of *Flavobacterium columnare*. *Fish Pathol.*, 34, 217-218.
- Welker, T. L., C. A. Shoemaker, C. R. Arias and P. H. Klesius (2005): Transmission and detection of *Flavobacterium columnare* in channel catfish *Ictalurus punctatus*. *Dis. Aquat. Org.*, 63, 129-138.
- Yeh, H. Y., C. A. Shoemaker and P. H. Klesius (2006): Sensitive and rapid detection of *Flavobacterium columnare* in channel catfish *Ictalurus punctatus* by a loop-mediated isothermal amplification method. *J. Appl. Microbiol.*, 100, 919-925.